

Joystick-controlled everyday mobility.

The JoyGlide is a lightweight powered wheelchair that is controlled with an easy-to-use joystick.

Completely customisable.

JoyGlide is a fully-scripted wheelchair, built to the occupant's precise requirements.

Built to take on whatever life throws at it.

JoyGlide is a versatile wheelchair designed to perform well in a variety of conditions, indoors and out. With a lightweight yet strongly engineered frame,

JoyGlide is highly manoeuvrable and easy to navigate through tight spaces like crowded cafés. The battery system and in-wheel motor are located on one side, making it easy to fold the chair for transport or storage.

Power in your hands.

With the simple movement of the joystick, the JoyGlide has the power to make travelling up hills and steep ramps effortless. The battery is powerful enough to take you on a distance of up to 15km on a single charge.

Designed by experience.

Today's JoyGlide is the result of decades of user feedback and product development to deliver the features that make JoyGlide a popular powered wheelchair.



Customised seating.

For ultimate comfort many seating options are available. Choose customised cushion, backrest and headrest.



Anti tip wheels.

For peace of mind, JoyGlide features swingaway anti tip wheels to avoid tipping on inclines.



Wheels.

The 16" or 24" drive wheels and 7" front castors have comfortable pneumatic tyres for an air-soft ride.



Electronics.

Fully programmable joystick control system.



Folding frame.

JoyGlide has the ability to be folded for transport.



Battery range.

The powerful battery gives you a range of up to 15km on a single charge.



Colour choices.

Choose from an extensive colour palette for the ultimate customisation.

The features that make JoyGlide such a joy.

Features.

- Easy-to-use joystick controls
- High quality in-wheel motors
- Low maintenance with accessible electronics
- Turning circle between walls 1145mm
- Max ramp incline/decline 10° (as per AS3695 testing)
- 16" or 24" drive wheel option with Quick Release
- Electromagnetic brakes and manual
- Full length standard armrests
- Swing-away footplates

Safety.

- Full compliance AS 3695
- Crash tested to AS 3696:19
- Max user weight 125kg

Arrange a test drive.

Chat to a Glide wheelchair specialist and experience the extraordinary features of the JoyGlide in person. Find your nearest distributor at glide.com.au

Specifications.



Recommended User Weight
Up to 125kg AS 3695 full compliance



Overall Length
1020mm



Seat Widths
37 to 52cm



Overall Width
600mm (wheel base)



Seat Depths
33 to 52cm



Seat Height
46cm (top of rail using 0-46 set up)



Back Height
1030mm Top of push handle (custom height available)



Arm Rests
Fixed



Wheelchair Weight
30.4kg



Leg Rests
Swing-away detachable



Frame Colour
Wide choice



Brakes
Electromagnetic and manual park



www.glide.com.au

(08) 9345 3400 | sales2@glide.com.au

25 Ledger Road, Balcatta, Western Australia 6021

@GlideProductsAustralia Glide Products